

MRI Breast Implant Routine Guideline (61A)

CPT Code: 77059

Indications: Implant pain, Implant Rupture

Contrast: No

| Sequence | Plane | Slice Thickness | Gap | FOV | Notes |
|-----------------|------------------|------------------------|------------|------------|---|
| T2 | Axial | 4mm | .6mm | 310-350 | Silicone Suppression, Use Silicone Peak |
| T2 | Axial | 4mm | .6mm | 310-350 | Water Suppression, Use Water Peak |
| TIRM | Sagittal Left | 4mm | .6mm | 310-350 | Water Suppression, Use Water Peak |
| TIRM | Sagittal Routine | 4mm | .6mm | 310-350 | Water Suppression, Use Water Peak |